

AGM Container Controls, Inc. 3526 E. Fort Lowell Rd., Tucson, AZ 85716

Tel: 520-881-2130

Email: sales@agmcontainer.com Website: www.agmcontainer.com

PRODUCT SPECIFICATIONS

AGM Part Number 920545

Molecular Sieve, Type 13, 8 x 12 Mesh Beaded

<u>CHEMICAL FORMULA</u> $Na_{86}[(AIO_2)_{86}(SiO_2)_{106}] \cdot nH_2O$

PRODUCT DESCRIPTION

Molecular sieve is a multiple purpose, highly porous, high capacity alkali metal alumino-silicate in the spherical form. It is the sodium form of the Type X crystal structure with pore diameters of approximately 10Å. It can adsorb all molecules that can be adsorbed by 3A, 4A, and 5A molecular sieve. Type 13X can also adsorb molecules, such as aromatics and branched-chain hydrocarbons, which have large critical diameters.

SPECIFICATIONS				
Property	Value		Unit	
Bead Size	8x12		Mesh	
Nominal Pore Opening	10		Angstroms	
Equilibrium Water Capacity @ 25°C	≥	26	%	
Heat of Adsorption	1800		BTU/lb of H ₂ O	
Bulk Density	≥	.64	g/ml	
	≥	40	lbs/ft ³	
Crush Strength	≥	8	lbs	
	≥	35	N	
Size Qualification	2	97	%	
Package Moisture	≤	1.5	% Weight	

TYPICAL APPLICATIONS

Removal of H_2O and CO_2 in air plant purification systems for air separation processes, as a desiccant in air dryers and PSA oxygen systems.

PACKAGING INFORMATION

Drums 308 lb (140 kg)

Super Sacks Fill weights to order specification

HANDLING AND STORAGE RECOMMENDATIONS

Store molecular sieve in a dry location to prevent premature water adsorption. Reseal packages after opening to prevent contamination and unintended water adsorption. We recommend that you rotate stock so oldest material is used first.

HEALTH & SAFETY INFORMATION

Health and Safety Information is available on the product MSDS.